

Brain Tumor Detection Using Deep Learning (ResNet-18 on MRI Images)

Code:

1):

%% Brain Tumor Classification using ResNet-18

% Author: Vijayagiri Shanthipriya

% Description: Trains a ResNet-18 model to classify brain tumors

% (glioma, meningioma, notumor, pituitary)

clc; clear; close all;

%% STEP 1: Set dataset folder path

**dataDir ='C:\Users\shanthipriya
vijayagiri\OneDrive\Desktop\Matlabproject\Braintumordetection';**

if ~isfolder(dataDir)

error('Data folder not found. Please check path.');

end

%% STEP 2: Create image datastore

**imds = imageDatastore('C:\Users\shanthipriya
vijayagiri \OneDrive\Desktop\Matlabproject\', ...
 'IncludeSubfolders', true, ...
 'LabelSource', 'foldernames');**

%% STEP 3: Display label distribution

tbl = countEachLabel(imds);

disp('Label counts:');

```
disp(tbl);
```

```
%% STEP 4: Check and fix grayscale images
```

```
inputSize = [224 224 3]; % ResNet-18 input size
```

```
imds.ReadFcn = @(filename) readAndFix(filename, inputSize);
```

```
function I = readAndFix(filename, inputSize)
```

```
    % Read image
```

```
    I = imread(filename);
```

```
    % Convert grayscale to RGB
```

```
    if size(I,3) == 1
```

```
        I = cat(3, I, I, I);
```

```
    end
```

```
    % Resize
```

```
    I = imresize(I, inputSize(1:2));
```

```
    % Convert to uint8 if needed
```

```
    if ~isa(I, 'uint8')
```

```
        I = im2uint8(I);
```

```
    end
```

```
end
```

```
%% STEP 5: Split dataset into training and validation
```

```
[imdsTrain, imdsVal] = splitEachLabel(imds, 0.8, 'randomized');
```

%% STEP 6: Load pretrained ResNet-18

```
net = resnet18;
```

%% STEP 7: Modify final layers for 4 classes

```
numClasses = numel(categories(imds.Labels));
```

```
lgraph = layerGraph(net);
```

```
newLayers = [
```

```
    fullyConnectedLayer(numClasses, 'Name','fc_new')
```

```
    softmaxLayer('Name','softmax_new')
```

```
    classificationLayer('Name','classoutput')];
```

```
lgraph = replaceLayer(lgraph,'fc1000',newLayers(1));
```

```
lgraph = replaceLayer(lgraph,'prob',newLayers(2));
```

```
lgraph =
```

```
replaceLayer(lgraph,'ClassificationLayer_predictions',newLayers(3));
```

%% STEP 8: Create augmented image datastores

```
augimdsTrain = augmentedImageDatastore(inputSize(1:2), imdsTrain,  
'ColorPreprocessing','gray2rgb');
```

```
augimdsVal = augmentedImageDatastore(inputSize(1:2), imdsVal,  
'ColorPreprocessing','gray2rgb');
```

%% STEP G: Define training options

```
options = trainingOptions('sgdm', ...
```

```
    'MiniBatchSize', 8, ...
```

```
'MaxEpochs', 8, ...  
'InitialLearnRate', 1e-4, ...  
'Shuffle', 'every-epoch', ...  
'ValidationData', augimdsVal, ...  
'ValidationFrequency', 30, ...  
'Verbose', true, ...  
'Plots', 'training-progress');
```

%% STEP 10: Train the network

```
trainedNet = trainNetwork(augimdsTrain, lgraph, options);
```

%% STEP 11: Save trained model

```
modelSaveFile = fullfile(pwd, 'trainedResNet18_brain.mat');  
save(modelSaveFile, 'trainedNet');
```

%% STEP 12: Evaluate performance

```
YPred = classify(trainedNet, augimdsVal);  
YValidation = imdsVal.Labels;
```

```
accuracy = mean(YPred == YValidation);  
fprintf('\nValidation accuracy: %.2f%%\n', accuracy * 100);
```

%% STEP 13: Confusion matrix

```
figure;  
plotconfusion(YValidation, YPred);
```

```
title('Confusion Matrix - Brain Tumor Classification');
```

```
2):
```

```
clc;
```

```
clear;
```

```
close all;
```

```
% Step 1: Load trained network
```

```
modelFile = 'trainedResNet18_brain.mat';
```

```
if isfile(modelFile)
```

```
    load(modelFile, 'trainedNet');
```

```
    disp('  Trained model loaded successfully.');
```

```
else
```

```
    error(' + Trained model file not found!');
```

```
end
```

```
% Step 2: Select image
```

```
[filename, pathname] = uigetfile({'*.jpg;*.png;*.jpeg'}, 'Select a Brain  
MRI Image');
```

```
if isequal(filename,0)
```

```
    disp('No image selected.');
```

```
    return;
```

```
end
```

```
imgPath = fullfile(pathname, filename);
```

```
img = imread(imgPath);
```

```
figure, imshow(img), title('Input Brain MRI Image');
```

% Step 3: Resize image to match network input

```
img_resized = imresize(img, [224 224]);
```

% Step 4: Predict using trained model

```
[predictedLabel, scores] = classify(trainedNet, img_resized);
```

```
disp(['📄 Predicted Tumor Type: ', char(predictedLabel)]);
```

% Step 5: Convert to grayscale for area detection

```
grayImg = rgb2gray(img);
```

% Step 6: Apply segmentation to highlight tumor region

```
bw = imbinarize(grayImg, 'adaptive');
```

```
bw = imfill(bw, 'holes');
```

```
bw = bwareaopen(bw, 80); % remove small noise
```

% Step 7: Calculate tumor area

```
tumorArea = bwarea(bw);
```

```
disp(['📄 Tumor Area (in pixels): ', num2str(tumorArea)]);
```

% Step 8: Display segmented tumor region

```
figure, imshow(img), title('Detected Tumor Region');
```

```
hold on;
```

```
visboundaries(bw, 'Color', 'r', 'LineWidth', 0.8);
```

```
hold off;
```

% Step G: Detect tumor presence clearly

thresholdArea = 500; % you can adjust this value

if tumorArea > thresholdArea

**msg = sprintf('🚨 Tumor Detected!\nType: %s\nArea: %.2f pixels',
 string(predictedLabel), tumorArea);**

else

msg = '✅ No Tumor Detected. Brain appears normal.';

end

disp(msg);

msgbox(msg, 'Tumor Detection Result');

disp('✅ Detection completed successfully.');